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- ・Coxsackieviruses A1 to A21, A24, and B1 to 6
- ・Echoviruses (enteric cytopathic human orphan viruses) 1 to 7, 9, 11 to 21, 24 to 27, and 29 to 33
- ・Enteroviruses 68 to 71, 73 to 91, and 100 to 101
- ・Polioviruses types 1 to 3

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DNA polymerases；AmpliTaq Gold® DNA Polymerase with Buffer II and MgCl₂

ゲル抽出；FastGene Gel/PCR Extraction Kit

Cycle sequence kits；BigDye® Terminator v3.1 Cycle Sequencing Kit

Dye Terminator Removal；FastGene Dye Terminator Removal Kit

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Briefly, 5 μ l of the RNA mixed with 1 μ l of 50 μ M of each primer (EVP-2, 5'-CCT CCG GCC CCT GAA TGC GGC TAA-3' and OL68-1, 5'GGT AAQ TTC CAC CAC CAY AA-3', Y is either C or T) was heated for 5 min at 95 °C and immediately put on ice. Then 43 μ l of reactive cocktail containing 5 μ l of 10X DNA polymerase reaction buffer, and 1 μ l of each of the following was added, 40 U of RNase inhibitor (Promega), 200 U MMLV (BRL), 2.5 U of cloned PfuTurbo DNA polymerase (Stratagene) and 15 mM of dNTP (Amershan). The one-step RT-PCR program consisted of incubation for 1 h at 37 °C, for 5 min at 94 °C, 40 cycles of incubation for 1 min at 95 °C, 30 sec at 55 °C, 1 min at 72 °C, and followed by 5 min at 72 °C. The second PCR was performed at a final volume of 50 μ l including 2.5 U of cloned PfuTurbo DNA polymerase, 1 μ l of 15 mM dNTP, 4.5 μ l 10X Pfu reaction buffer, and 1 μ l of 50 μ M of each primer (EVP-4, 5'-CTA CTT TGG GTG TCC GTG TT-3' and OL68-1) in the same PCR program. A fragment of approximately 650 base pairs spanning the 5'-noncoding region to one third of VP2 including the entire region of VP4 was amplified. Cycle sequencing was performed using the purified PCR products with the Prism Ready Reaction Dideoxy Terminator cycle sequencing kit (Perkin-Elmer Corporation-Applied Biosystems, Foster City, CA, USA) and the 207bp nucleotide sequence was determined with an automated sequencer (ABI Model 373A).

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Cycle Sequencing on the System 9700, 9600, 2700, or 2400

To sequence single- and double-stranded DNA on the GeneAmp® PCR System 9700 (in 9600 emulation mode), 9600, or 2400:

1.	Place the tubes in a thermal cycler and set to the correct volume.
2.	Perform an initial denaturation. <ol style="list-style-type: none"> a. Rapid thermal ramp to 96 °C b. 96 °C for 1 min
3.	Repeat the following for 25 cycles: <ul style="list-style-type: none"> • Rapid thermal ramp to 96 °C • 96 °C for 10 sec • Rapid thermal ramp to 50 °C • 50 °C for 5 sec • Rapid thermal ramp to 60 °C • 60 °C for 4 min
4.	Rapid thermal ramp to 4 °C and hold until ready to purify.
5.	Spin down the contents of the tubes in a microcentrifuge.
6.	Proceed to Chapter 4, "Purifying Extension Products."

*Rapid thermal ramp is 1 °C/second.

Table3 Cycle sequence standard protocol

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Performing Cycle Sequencing To perform cycle sequencing on BAC DNA:

1.	Place the tubes in a thermal cycler and set the volume to 20 μ L.
2.	Heat the tubes at 95 $^{\circ}$ C for 5 minutes.
3.	Repeat the following for 30 cycles: [*] <ul style="list-style-type: none"> • Rapid thermal ramp[†] to 95 $^{\circ}$C • 95 $^{\circ}$C for 30 seconds. • Rapid thermal ramp to 50–55 $^{\circ}$C (depending on template) • 50–55 $^{\circ}$C for 10 seconds. • Rapid thermal ramp to 60 $^{\circ}$C • 60 $^{\circ}$C for 4 minutes.
4.	Rapid thermal ramp to 4 $^{\circ}$ C and hold until ready to purify.
5.	Spin down the contents of the tubes in a microcentrifuge.
6.	Proceed to Chapter 4, "Purifying Extension Products."

^{*}Some laboratories have found that increasing the number of cycles gives better results.

[†]Rapid thermal ramp is 1 $^{\circ}$ C/sec.

Table4 Cycle sequence BAC DNA protocol

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DNA Quantity

Quantitating DNA If possible, quantitate the amount of purified DNA by measuring the absorbance at 260 nm or by some other method.

Template Quantity

The table below shows the amount of template to use in a cycle sequencing reaction.

Template	Quantity
PCR product:	
100–200 bp	1–3 ng
200–500 bp	3–10 ng
500–1000 bp	5–20 ng
1000–2000 bp	10–40 ng
>2000 bp	20–50 ng
Single-stranded	25–50 ng
Double-stranded	150–300 ng
Cosmid, BAC	0.5–1.0 μ g
Bacterial genomic DNA	2–3 μ g

Note: In general, higher DNA quantities give higher signal intensities.

The template quantities stated above should work with all primers. You may be able to use even less DNA when using capillary instruments for detection. The amount of PCR product to use in sequencing also depends on the length and purity of the PCR product.

Using BigDye Terminator v1.1/3.1 Sequencing Buffer

The BigDye[®] Terminator v1.1/3.1 Sequencing Buffer (5X)^{*} is supplied at a 5X concentration. If you use it for sequencing reactions, be sure the final reaction volume is at a concentration of 1X. For example, for a half reaction in 20 μ L final volume, you would use 4 μ L of ready reaction premix and 2 μ L of BigDye sequencing buffer as shown below.

Reagent	Concentration	Volume
Ready Reaction Premix	2.5X	4 μ L
BigDye Sequencing Buffer	5X	2 μ L
Primer	–	3.2 pmol
Template	–	See "Template Quantity" on page 2-6.
Water	–	to 20 μ L
Final Volume	1X	20 μL

Note: The use of this buffer without optimization may result in deterioration of sequence quality. Applied Biosystems does not support diluted reactions or guarantee the performance of BigDye[®] chemistry when it is diluted.

Table5 Cycle sequence DNA Quantity

Table6 Reaction mixtures

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